



## Enterprise Low Bed-NEW!!


### MMC has new low beds Enterprise 9000x from ARJO

**Beginning May 18<sup>th</sup>, 2021** all low beds ordered will be this new bed



#### What you need to know:

- Bed exit alarm cords connect to our call light system just like our regular beds
  - These special connection cables come with the bed
  - Contact MOWS if cable is not provided (8-3937 or page 6115). After hours contact House Supervisor.
- Must zero bed when patient is out of the bed and then weigh the patient for the sensitivity levels of the alarm to work correctly
- To set bed alarm, must hold down the *In bed* or *Egress* button for the sensitivity level desired.



- There are 10 alarm sensitivity levels
- Weight capacity is 550 pounds

The ordering process will not change. Continue to place a *Specialty Bed Referral* and specify need for low bed.

**See Quick Reference Guide on the Enterprise controls posted on Memorial Central. Additional resource video located on Workplace.  
Contact your Nurse Manager with any questions.**